


BioBiteCorrector Sleep Appliance

The anti-snoring device with a titanium ball-joint hinge



starts at € 449,- *

Advantages

- Anti-snoring hinge made completely of titanium
- Biocompatible and hypoallergenic
- Low-level vertical bite opening
- Direct adjustability
- High wearing comfort thanks to ball joints
- Digital 3D design
- Splint surface design with occlusal support
- Reliable BioBiteCorrector hinge benefits
- Premium-quality milled splints
- High wearing comfort thanks to flat design
- High effectiveness

BBC-Orthotec GmbH
Innsbrucker Str. 2
83435 Bad Reichenhall

Tel 08651 9650099
Fax 08651 9650098
E-Mail: lab@bbc-orthotec.de

* All prices excl. 7% VAT,
valid from 1 Jan. 2019.

BioBiteCorrector® Herbst

The custom-cast splint appliance



starts at € 449,- *

BBC-Orthotec LAB is a laboratory specialising in devices for orthodontics and dental sleep medicine. We manufacture orthodontic and anti-snoring appliances based on our innovative products.

Digital workflow:

We manufacture our devices in Germany, using intraoral scans or plaster models and both digital and conventional techniques.

The **BioBiteCorrector Sleep Appliance** is designed digitally, using 3D software. Then the splints are milled (CAD/CAM process).

The **BioBiteCorrector Herbst Appliance** is designed using 3D software (CAD). In a 3D printing process, the bands including the BBC attachments are printed in castable resin and then cast in metal.

BioBiteCorrector Herbst Appliance

The custom-cast Herbst appliance



starts at € 449,- *

BioBiteCorrector Twin appliance

The modern functional appliance



starts at € 379,- *



Advantages

- Digital 3D design
- High wearing comfort thanks to flat BBC SA design
- Construction bite is not required
- All BBC components can be replaced chairside
- Reliable BioBiteCorrector hinge benefits

Advantages

- **The BioBiteCorrector MS & SP can be unscrewed from the Schwartz plate and later fixed to a multi-bracket appliance on the same patient**
- Construction bite is not required
- Without lateral occlusal pads, for easy deep bite correction and higher wearing comfort
- Does not cause any tongue space restriction, so that speaking comfort is higher
- Also available as a milled splint
- Reliable BioBiteCorrector hinge benefits