


BioBiteCorrector® Sleep Appliance

The anti-snoring device with
a titanium ball-joint hinge



Advantages:

- Anti-snoring hinge made completely of titanium
- Biocompatible and hypoallergenic
- Thermoactive flexibility
- MMA-free (0% residual monomer)
- Self-adjusting with memory effect
- Low-level vertical bite opening
- Infinitely variable calibration
- Calibrated activation with spacer
- High wearing comfort thanks to ball joints
- Digital 3D-Design (CAD/CAM)
- Splint surface design with occlusal support
- Reliable BioBiteCorrector® hinge benefits
- Premium-quality milled splints
- High wearing comfort thanks to flat design

Content controlled, updates and corrections are subject to change without prior notice - Bild Titelseite © shutterstock.com - Edited 120 BBC LAB EN 04/2020

orthotec
bbc LAB
Orthodontic & Dental Sleep Medicine Laboratory

BBC-Orthotec GmbH
Innsbrucker Str. 2
D-83435 Bad Reichenhall

Phone: +49 8651 9650099
Fax: +49 8651 9650098
E-Mail: lab@bbc-orthotec.de

orthotec
bbc LAB

Your orthodontic
and dental sleep
laboratory



*All prices excl. VAT, valid from 1st Jan. 2020.

BBC-Orthotec LAB is a laboratory specialising in devices for orthodontics and dental sleep medicine. We manufacture orthodontic and anti-snoring appliances based on our innovative products.

Digital workflow:

We manufacture our devices in Germany, using intraoral scans or plaster models and both digital and conventional techniques.

The BioBiteCorrector® **Sleep Appliance** is designed digitally, using 3D software. Then the splints are milled (CAD/CAM process).

The BioBiteCorrector® **Herbst Appliance** is designed digitally, using 3D software. In a Metal-3D printing process (CAM / laser melting) the appliance is printed.



BioBiteCorrector® Herbst Appliance

The CAD/CAM BBC Herbst



starts at **499,- €***

Advantages:

- Digital 3D design
- Metal-3D printing
- High wearing comfort thanks to flat BBC SA design
- Construction bite is not required
- All BBC components can be replaced chairside
- Reliable BioBiteCorrector® hinge benefits

BioBiteCorrector® Twin Appliance

The modern functional appliance



starts at **359,- €***

Advantages:

- The BioBiteCorrector® MS & SP can be unscrewed from the Schwartz plate and later fixed to a multi-bracket appliance on the same patient
- Construction bite is not required
- Without lateral occlusal pads, for easy deep bite correction and higher wearing comfort
- Does not cause any tongue space restriction, so that speaking comfort is higher
- Also available as a milled splint
- Reliable BioBiteCorrector® hinge benefits