



BBC-Orthotec GmbH
Innsbrucker Str. 2
D-83435 Bad Reichenhall

Phone: +49 8651 9650099
Fax: +49 8651 9650098
E-Mail: lab@bbc-orthotec.de



BioBiteCorrector[®] Herbst

The fixed
dental brace for
overbite correction




BioBiteCorrector® Herbst

The fixed dental brace for overbite correction



The BioBiteCorrector® Herbst is a dental brace, also known as a functional orthodontic appliance, and the modern version of the well-known Herbst hinge.

This version of the BioBiteCorrector® Herbst is a sturdy and fixed appliance designed to correct an overbite. An overbite means your lower jaw is too far back, relative to your upper jaw, and your teeth do not bite properly.



improper bite

proper bite

Advantages:

- Quick Result
- The appliance does not require any active cooperation and works 24 hours a day; this is why it is so successful.
- Modern digital CAD/CAM technology
- The extremely flat and slender design ensures high wearing and speaking comfort. The metal base makes the appliance very sturdy.
- Biocompatible titanium
- The BBC telescope is made entirely of titanium and laser-welded.
- Proper bite
- The BBC moves your lower jaw in the right position, so that your teeth bite properly.
- Shorter bracket wearing time
- The appliance corrects your lower jaw position before treatment with a brace fixed by brackets. This will shorten the treatment time with brackets.

